



Recombinant *Saccharomyces cerevisiae* Protein YRO2 (YRO2)

Product Code	CSB-CF339911SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P38079
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MSDYVELLKRGGNEAIKINPPTGADFHITSRGSDWLFTVFCVNULLFGVILVPLM FRKPVK DRFVYYTAIAPNLFMSIAYFTMASNLGWIPVRAKYNHVQTSTQKEHPGYRQIFY ARYVGW FLAFPWPIIQMSLLGGTPLWQIAFNVGMTEIFTVCWLIAACVHSTYKWGYTIGI GAAIV VCISLMTTTFNLVKARGKDVSNVFITFMSVIMFLWLIAYPTCFGITDGGNVLQPD SATIF YGIIDLLILSILPVLFMPLANYLGIERLGLIFDEEPAEHVGPVAEKKMPSPASFKSS DSD SSIKEKLLKLLKHHKDKKKAKKAKKAKKAKKAQEEEEEDVATDSE
Source	in vitro E.coli expression system
Target Names	YRO2
Protein Names	Recommended name: Protein YRO2
Expression Region	1-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.