



Recombinant Pan troglodytes Trace amine-associated receptor 6 (TAAR6)

Product Code	CSB-CF719640EQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5W8W0
Product Type	Transmembrane Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Sequence	MSSNSSLLVAVQLCYPNVNGSCVETLYSPGSRVILYIVFGFGAVLAVFGNLLVM ISILHF KQLHSPTNFLVASLACADFLVGVTVMPPFSMVRTVESCWYFGRSFCTFHTCCD VAFCYSSL FHLCFISIDRYIAVTDPLVYPTKFTVSVSGICISVSWILPLMYSGAVFYTG VYDDG LEEL SDALNCIGGCQTVVNQNWVLIDCLSFFIPTFIMIILYGNIFLVARRQAKKIENTGS KTES SSESYKARVARRERKAAKTLGVTVVAFMISWLPYSIDSLIDAFMGFITPAYIYEIC CWCA YYNSAMNPLIYALFYPWFRKAIKVIVTGQVLKNSSATMNLFSEHI
Source	in vitro E.coli expression system
Target Names	TAAR6
Protein Names	Recommended name: Trace amine-associated receptor 6 Short name= TaR-6 Short name= Trace amine receptor 6 Alternative name(s): Trace amine receptor 4 Short name= TaR-4
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.