



Recombinant Mouse Probable G-protein coupled receptor 139 (Gpr139)

Product Code	CSB-CF800329MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q80UC8
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MEHTHAHLAANSSACGLGFVPVYYSFLLCLGLPANILTVIILSQLVARRQKSSY NYLLA LAAADILVLFVIVDFLLEDFILTMQMPLIPDKIIEVLEFSSIHTSIWITVPLTVDRY AVCHPLKYHTVSYPARTRKVILSVYITCFLTSSIPYYWWPNIWTEYISTSMHHVL VWIHC FTVYLVPCSIFFILNSIIVYKLRKSNFRLRGYSTGKTTAILFTITSIFATLWAPRIIM LYHLYGAPIQNPWL VHIMLDVANMLALLNTAINFFLYCFISKRFRTMAAATLKAL FKCQK QPVQFYTNHNFSITSSPWISPANSHCIKMLVYQYDKHKGPIKVSP
Source	in vitro E.coli expression system
Target Names	Gpr139
Protein Names	Recommended name: Probable G-protein coupled receptor 139 Alternative name(s): G(q)-coupled orphan receptor GPRg1 G-protein-coupled receptor PGR3
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.