



Recombinant Human Olfactory receptor 2K2 (OR2K2)

Product Code	CSB-CF844071HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGT1
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNITLYSFLYVLRLNEGNSTGRNIDERKKMQGENFTIWSIFFLEGFSQYPGLEV VLFVFS LVMYLTTLGNSTLILITILDSRLKTPMYLFLGNLSFMDICYTSASVPTLLVNLLSS QKT IIFSGCAVQMYLSLAMGSTECVLLAVMAYDRYVAICNPLRYSIIMNRCVCARMA TVSWVT GCLTALLETSFALQIPLCGNLIDHFTCEILAVLKLACTSSLLMNTIMLVVSILLPI ML LVCISYIFILSTILRITSAEGRNKAFTSCGAHLTVVILYYGAALSMYLPSSSNAQK IDK IISLLYGVLTPLNPIIYSLRNKEVKDAMKLLGKITLHQTHEHL
Source	in vitro E.coli expression system
Target Names	OR2K2
Protein Names	Recommended name: Olfactory receptor 2K2 Alternative name(s): HTPCRH06 Olfactory receptor OR9-17
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.