



Recombinant Serpentine receptor class T-55 (srt-55)

Product Code	CSB-CF327342CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34571
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	ICFDLKLQCWPMEIQEMALMLTARNTARYDCSGKSKSEWYETGQKRLGWGI YYISSGLF FQLIGWPVIWVFITKFSMTNALKVYRIMVFIGLIEITEIWGNSVFPGFVAVFGEVY CTSP ILMTIVGKMTMVQWVLGSSSAFLGFHRLCDMIQKLEWLVTNTKTGLWLTVL FFYACYG SIFDFTVLFNSDYMAPLLDPMIGKQGIYSNNFLYFHNIIVATTLILVYACLCTLWS SRE MNTSSLHVSKFQRSILLQSICISLTYAIPAISFVTMFVLPKWFHVSVDITYQLSG GLP FIMYICLNKRVREEFLHLLRVCRKAEKSQVAVIPLGNSTVSAFNN
Source	in vitro E.coli expression system
Target Names	srt-55
Protein Names	Recommended name: Serpentine receptor class T-55 Short name= Protein srt-55
Expression Region	19-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.