



Recombinant Chicken Melatonin receptor type 1C

Product Code	CSB-CF344189CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49288
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	MERPGSNGSCSGCRLEGGPAARAASGLAAVLIVTIVVDVLGNALVILSVLRNKK LRNAGN IFVVLSVADLVVAVYPYPLILSAIFHNGWTMGNIHCQISGFLMGLSVIGSIFNITA IAI NRYCYICHSLRYDKLFNLKNTCCYICLTWTLTVVAIVPNFFVGSQYDPRIYSCT FAQTV STSYTITVVVHFIVPLSIVTFCYLRIWILVIQVKHRVRQDCKQKIRAADIRNFLTM FVV FVLFAVCWGPLNFIGLAVSINPSKVQPHIPEWLFVLSYFMAYFNSCLNAVIYGLL NQNFR KEYKRILLMLRTPRLLFIDVSKGGTEGLKSKPSPAVTNNNQAEIHL
Source	in vitro E.coli expression system
Protein Names	Recommended name: Melatonin receptor type 1C Short name= Mel-1C-R Short name= Mel1c receptor Alternative name(s): CKB
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.