



Recombinant Gorilla gorilla gorilla fMet-Leu-Phe receptor (FPR1)

Product Code	CSB-CF008854GGZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79176
Product Type	Transmembrane Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Sequence	NSSLPTNISGGTPAVSAGYLFLDIVTYLVFAVTFVLGVLGNGLVIWVAGFRMTH TVTTIS YLNLA VADFCFTSTL PFFMVKKAMGGHWPFGWFLCKFIFTIVDINLFGSVFLIAL IALDR CVCVLHPVWTQNHRTVSLAKKVIIGPWVMALLLTLPVIIRVTTVPGKMGTVACT FNFSPW TNDPKERIKVAVAMLTVRGIIRFIIGFSAPMSIVAVSYGLIATKIDKQGLIKSSRTL RVL SFVAAAFFLSWSPYQVVALIATVRIRELLQGMYKEIGIAVDVTSALAFFNSCLNP MLYVF MGQDFRERLIHALPASLERALTEDSTQTS DTATNSTLPSAEVALQA
Source	in vitro E.coli expression system
Target Names	FPR1
Protein Names	Recommended name: fMet-Leu-Phe receptor Short name= fMLP receptor Alternative name(s): N-formyl peptide receptor Short name= FPR N- formylpeptide chemoattractant receptor
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.