



Recombinant Danio rerio Membrane progesterin receptor gamma-B (paqr5b)

Product Code	CSB-CF728495DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6DC77
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MLSLIKLQRVFNHQPFAFHEDGIISGYRHPRSSATECVWSLFQLTNETLNV WTHFLPT WYFLWKLMTVLLMEDVWNEAYTWPLLVFLFSCCVYPLASSCAHTFSSMSTRA RHICYFFD YGALSFYSLGSAISYSAYVFPDAWLSSSFHAYYISVAVFNTVLSTSLACYSRLG LPLLHY SHDIVERFSERQCPRMSKVLRLAFAYPYLFDNIPLFYRLFVFCVGEVCTDNEAN SVHVQH TLLAFLTSFLFATHLPERLAPGRFDYIGHSHQLFHVCAIIGTHFQMKAEMDMGL RRSQL LASAPAISFNNTIGAALLCVSVSLGIICVYSLPLLYSSNPKN TANKE
Source	in vitro E.coli expression system
Target Names	paqr5b
Protein Names	Recommended name: Membrane progesterin receptor gamma-B Short name= mPR gamma-B Alternative name(s): Progesterin and adipoQ receptor family member V-B
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.