



Recombinant Human Olfactory receptor 1L6 (OR1L6)

Product Code	CSB-CF836760HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGR2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSYFYRLKLMKEAVLVKLPFTSLPLLLQTLRKSIRDMEIKNYSSSTSGFILLGLS SNPQL QKPLFAIFLIMYLLAAVGNVLIIPAIYSDPRLHTPMYFFLSNLSFMDICFTTVIVPK MLV NFLSETKVISYVGCLAQMYFFMAFGNTDSYLLASMAIDRLVAICNPLHYDVVMK PRHCLL MLLGSCSISHLHSLFRVLLMSRLSFCASHIIKHFFCDTQPVLKLSCSDTSSSQM VVMETET LAVIVTPFLCIIFSYLRLIMVTVLRIPSAAGKWKAFSTCGSHLTAVALFYGSIIVVYF RPL SMYSVVRDRVATVMYTVVTPMLNPFYISLRNKDMKRGLKKLQDRIYR
Source	in vitro E.coli expression system
Target Names	OR1L6
Protein Names	Recommended name: Olfactory receptor 1L6 Alternative name(s): Olfactory receptor 1L7 Olfactory receptor OR9-30
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.