



Recombinant Human Olfactory receptor 2M2 (OR2M2)

| | |
|--------------------------|--|
| Product Code | CSB-CF850420HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q96R28 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Homo sapiens (Human) |
| Sequence | MAWENQTFNSDFILLGIFNHSPHTFLFFLVLGIFLVAFMGNSVMVLLIYLDLDTQL HTPMY FLLSQLSLMDLMLICTTVPKMAFNLYSGSKSISMAGCVTQIFFYISLSGSECFLL AVMAY DRYIAICHPLRYTNLMNPKICGLMATFSWILGSTDGIIDAVATFSFSFCGSREIAH FFCE FPSLLILSCNDTSIFEEVIFICCIVMLVFPVAIIIASYARVILAVIHMGSGEGRCKAF TT CSSHLMVVGMYGGAALFMYIRPTSDHSPTQDKMVSVFYILTTPMLNPLIYSLRN KEVTRA FMKILGKKGKSESELPHKLYVLLFAKFFFLISIFFYDVKILALIMYIA |
| Source | in vitro E.coli expression system |
| Target Names | OR2M2 |
| Protein Names | Recommended name: Olfactory receptor 2M2 Alternative name(s): OST423 |
| Expression Region | 1-347 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |