



Recombinant Human Olfactory receptor 13C3 (OR13C3)

Product Code	CSB-CF850900HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGS6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MIVQLICTVCFLAVNTFHVRSSFDLKAADDMGEGINQTLVSEFLLGLSGYPKIEIV YFAL ILVMYLVILIGNGVLIASIFDSHFHTPMYFFLGNL SFLDICYTSSSV PSTLVSLISK KR NISFSGCAVQMFFGFAMGSTECLLLGMMAFDRYVAICNPLRYPILSKVAYVLM ASVSWL SGGINS AVQTLLAMRLPFCGNNIINHFACEILAVLKLACADISLNIITMVISNMAFL VLP LMVIFFSYMFILYTI LQMNSATGRRKAFSTCSAHLTVIIFYGTIFFMYAKPKSQD LIGE EKLQALDKLISLFYGVVTPMLNPILYSLRNKDVKA AVKYLLNKKPIH
Source	in vitro E.coli expression system
Target Names	OR13C3
Protein Names	Recommended name: Olfactory receptor 13C3 Alternative name(s): Olfactory receptor OR9-8
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.