



Recombinant *Loxodonta africana* Rhodopsin (RHO)

Product Code	CSB-CF019681LJC
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q68J47
Product Type	Transmembrane Protein
Immunogen Species	<i>Loxodonta africana</i> (African elephant)
Sequence	MNGTEGPNFYVPFSNKTGVVRSFPFEYPQYYLAEPWQFSMLAAYMFLLVLGF PINFLTLY VTVQHKNVRTPLNYILLNLAVANHFVFGGFTTTLTSLHGYFVFGSTGCNLE GFFATLG GEIALWSLVVLAIERVYVVCKPMSNFRFGENHAIMGVAFTWVMALACAAPPLV GWSRYIP EGMQCSCGIDYYTLKPEVNNESFVIYMFVVHFTIPMTIIFFCYGQLVFTVKEAAA QQQES ATTQKAEKEVTRMVIIMVIAFLICWVPYASVAFYIFTHQGSDFGPILMTLPAFFAK SSAI YNPVIYIMMNKQFRNCMLTTICCGKNPFGEEEGSTTASKTETSQVAPA
Source	in vitro E.coli expression system
Target Names	RHO
Protein Names	Recommended name: Rhodopsin
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.