



Recombinant Human Olfactory receptor 2T4 (OR2T4)

Product Code	CSB-CF854148HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NH00
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDNITWMASHTGWSDFILMGLFRQSKHPMANITWMANHTGWSDFILLGLFRQ SKHPALLC VVIFVVFLMALSGNAVLILLIHCD AHLHTPMYFFISQLSLMDMAYISVTVPKMLLD QVMG VNKISAPECGMQMFFYVTLAGSEFFLLATMAYDRYVAICHPLRYPVLMNHRVC LFLSSGC WFLGSVDGFTFTPITMTFFPFRGSREIHHFFCEVPAVLNLSCSDTSLYEIFMYLC CVLMLL IPVVISSSYLLILLTIHGMNSAEGRKKAFATCSSHLTVVILFYGAAIYTYMLPSSY HTP EKDMMVSVFYTILTPVVNPLIYSLRNKDVMGALKKMLTVEPAFQKAME
Source	in vitro E.coli expression system
Target Names	OR2T4
Protein Names	Recommended name: Olfactory receptor 2T4 Alternative name(s): Olfactory receptor OR1-60
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.