



Recombinant Phoca vitulina Rhodopsin (RHO)

Product Code	CSB-CF019681EUR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O62794
Product Type	Transmembrane Protein
Immunogen Species	Phoca vitulina (Harbor seal)
Sequence	MNGTEGPNFYVPFSNKTGVVRSPPFEFPQYYLAEPWQFSMLAAYMFLIVLGFP INFLTLY VTVQHKKLRTPLNILLNLAVADLFMVFGGFTTTLYTSLHGYFVFGPTGCNLEG FFATLG GEIALWLSLVLAIERVYVVCKPMSNFRFGENHAIMGVGFTWVMALACAAPPLV GWSRYIP EGMQCSCGIDYYTLKPEVNNESFVIYMFVVHFTIPMIVIFFCYGQLVFTVKEAAA QQQES ATTQKAEKEVTRMVIIMVIAFLICWVPYASVAFYIFTHQGSNFGPIFMTLPAFFAK AASI YNPVIYIMMNKQFRTCMITTLCCGKNPLGDDEVSAASKTETSQVAPA
Source	in vitro E.coli expression system
Target Names	RHO
Protein Names	Recommended name: Rhodopsin
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.