



Recombinant Rat Rhodopsin (Rho)

Product Code	CSB-CF019681RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51489
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MNGTEGPNFYVPFSNITGVVRSQFEQYLAEPWQFSMLAAYMFLIVLGF INFLTY VTVQHKKLRTPVLLNLAVALDFMVFVGGFTTLYTSLHGYFVFGPTGCNLEG FFATLG GEIGLWVSLVLAIRYVVVCKPMSNFRFGENHAIMGVAFTWVMALACAAPPLV GWSRYIP EGMQCSCGIDYYTLKPEVNNESFVIYMFVHFTIPMIVIFFCYGQLVFTVKEAAA QQQES ATTQKAEKEVTRMVIIMVIFFLICWLPYASVAMYIFTHQGSNFGPIFMTLPAFFA KTASI YNPIIYIMMNKQFRNCMLTSLCCGKNPLGDDEASATASKTETSQVAPA
Source	in vitro E.coli expression system
Target Names	Rho
Protein Names	Recommended name: Rhodopsin
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.