



Recombinant Macaca mulatta (Rhesus macaque) N-formyl peptide receptor 3 (FPR3)

Product Code	CSB-CF008856MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79191
Product Type	Transmembrane Protein
Immunogen Species	Macaca mulatta (Rhesus macaque)
Sequence	METNFSIPLNETEEVLPEPAGHTVLWIFSLLVHGVTFIFGVLGNGLVIWVAGFR MTRTVN TICYLNLALADFSFSAILPFHMOVSVAMREKWPFGTFLCKLVHVMIDINLFVSVYLI TIIA LDRICVLHPAWAQNHRTMSLAKRVMTGLWILTIVLTLPNFIFWTTISTTNGDTY CIFNY PFWGDTVVERMNVFITMAKVSLILHFIIGFSIPMSIITVCYGIIVAKIHKRMTKSS RPL HIFTAVVASFFICWFPYELTGILMAVWLKEILLNGKYKIILVLINPTSSLAFFNSCL NPS LYVFMGHNFQERLIRSLPTSLERALTEVPDSAQTSNTHHTTSASPPEETE
Source	in vitro E.coli expression system
Target Names	FPR3
Protein Names	Recommended name: N-formyl peptide receptor 3 Alternative name(s): FMLP-related receptor II Short name= FMLP-R-II Formyl peptide receptor-like 2
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.