



Recombinant Human Olfactory receptor 52I2 (OR52I2)

Product Code	CSB-CF844079HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NH67
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MCQQILRDCILLIHHLCINRKKVSLVMLGPAYNHTMETPASFLLVGIPGLQSSHL WLAIS LSAMYIALLGNTIIVTAIWMDSTRHEPMYCFLCVLAAVDIVMASSVVPKMVSIFC SGDS SISFSACFTQMFFVHLATAVETGLLLTMAFDRYVAICKPLHYKRILTPQVMLGM SMAITI RAIIAITPLSWMVSHLPFCGSNVVVHSYCEHIALARLACADVPSSLYSLIGSSL MVGSD VAFIAASYILILKAVFGLSSKTAQLKALSTCGSHVGMALYYLPGMASIYAAWLG QDVVP LHTQVLLADLYVIIPATLNPIIYGMRTKQLRERIWSYLMHVLFDHNSLGS
Source	in vitro E.coli expression system
Target Names	OR52I2
Protein Names	Recommended name: Olfactory receptor 52I2 Alternative name(s): Olfactory receptor OR11-12
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.