



Recombinant Dog B1 bradykinin receptor (BDKRB1)

Product Code	CSB-CF857721DO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BDQ5
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sequence	MASRAPLELLPLNRSQSPNATTCDDAPEAWDLLHRVLPSEVIIIICVCGLLGNL LVLAV LLRPRRRLNVAEMYLANLAASDLVFLGLPFWAANISNQFRWPFGLLCRLVN GVIKANL FISIFLVVAISRDRYRALVHPMATRRRRQARATCVLIWVAGSLLSVPTFLFRSIE AVPEL NNDACVLLHPPGAWHVARMVELNVLGFLPLAAIVFFNCHILASLRGRPEVR GARCGGP PDGRTTALILTFVAAFLVCWTPYHFFAFLEFLTQVQVVRGCFWENFKDLGLQY ASFFAFI NSCLNPVIYVFGRLFRTRVWDLFKQCAPRRPPAVSWSHRKRVLQLFWQN
Source	in vitro E.coli expression system
Target Names	BDKRB1
Protein Names	Recommended name: B1 bradykinin receptor Short name= B1R Short name= BK-1 receptor
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.