



Recombinant Human Putative transmembrane protein 185B (TMEM185B)

Product Code	CSB-CF887995HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9H7F4
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNPRGLFQDFNPSKFLIYTCLLLFSVLLPLRLDGIIQWSYWAVFAPIWLWKLLVV AGASV GAGVWARNPRYRTEGEACVEFKAMLIAVGIHLLLLMFEVLVCDRVERGTHFWL LVFMPLF FVSPVSVAACVWGFRHDSLELEILCSVNILQFIFIALKLDRIIHWPLVVFVPL WILMS FLCLVVLYYIVWSLLFLRSLDVVAEQRRTHVTMAISWITIVPLLTFEVLLVHRLD GHNT FSYVSIFVPLWLSLLTLMATTFRRKGGNHHWWFGIRRDFCQFLLEIFPFLREYGN ISYDLH HEDESDAEETSVPEAPKIAPIFGKKARVVITQSPGKYVPPPKLNIDMPD
Source	in vitro E.coli expression system
Target Names	TMEM185B
Protein Names	Recommended name: Putative transmembrane protein 185B Alternative name(s): Protein FAM11B
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.