



Recombinant Mouse C-X-C chemokine receptor type 1 (Cxcr1)

Product Code	CSB-CF800386MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q810W6
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MAEAEYFIWTNPEGDFEKEFGNITGMLPTGDYFIPCKRVPITNRQALVVFYALV SLLSLL GNSLVMLVILYRRRTRSVM DVYVLNLA IADLLFSLTLPFLAVSKLKGWIFGTPLC KMSVSL LKEFNFFSGILLACISVDRYLAI VHATR TLARKRYLVKFVCGI WGLSLILSLPFA IFR QAYKPFRSGTVCYEVLGEATDFRMTLRGLSHIFGFLPLLTMLVCYGLTLRML FKTHMR QKHRAMGVIFAVVLVFLCCLPYNLVLLSDTLLGAHLIEDTCERRNDIDQALYIT EILGF SHSCLNPIIYAFVGQNF RHEFLKILANHGLVRKEVLTHRRVAFHTSLTAIY
Source	in vitro E.coli expression system
Target Names	Cxcr1
Protein Names	Recommended name: C-X-C chemokine receptor type 1 Short name= CXC-R1 Short name= CXCR-1 Alternative name(s): High affinity interleukin-8 receptor A Short name= IL-8R A CD_antigen= CD181
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.