



Recombinant Mouse Leukotriene B4 receptor 1 (Ltb4r)

Product Code	CSB-CF013221MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O88855
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MAANTTSPAAPSSPGGMSLSLLPIVLLSVALAVGLPGNSFVVWSILKRMQKRT VTALLVL NLALADLAVLLTAPFFLHFLARGTWSFREMGCRLLCHYVCGISMYASVLLITIMSL DRSLA VARPFMSQKVRTKAFARWVLAGIWVVSFLLAIPVLVYRTVKWNNRTLICAPNY PNKEHKV FHLLFEAITGFLLPFLAVVASYSIDIGRRLQARRFRRSRRTGRLVVLILAFAAF PYHL VNLVEAGRTVAGWDKNSPAGQRLRLARYVLIALAFLLSSVNPVLYACAGGGLL RSAGVGF VVKLLEGTGSEVSSTRGGTLVQTPKDTPACPEPGPTDSFMTSSTIPESSK
Source	in vitro E.coli expression system
Target Names	Ltb4r
Protein Names	Recommended name: Leukotriene B4 receptor 1 Short name= LTB4-R 1 Short name= LTB4-R1
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.