



Recombinant Sargocentron diadema Rhodopsin (rho)

Product Code	CSB-CF019681SCO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79898
Product Type	Transmembrane Protein
Immunogen Species	Sargocentron diadema (Crown squirrelfish) (Adioryx diadema)
Sequence	MNGTEGPFYIPMVNTTGIVRSPYEYPQYYLVNPAAYAILGAYMFFLIIVGFPVN FMTLY VTLEHKKLRTPLNYILLNLAVADLFMVIGGFTTTTMYTSMHGYFVLGRLGCNLEG FFATLG GMISLWVSLAVLAIERWVVVCKPISNFRFGENHAIMGVSLTWGMALACTVPPLV GWSRYIP EGMQCSCGIDYYTRAEGFNNEFVLYMFFCHFTIPLTIIFFCYGRLLCAVKEAA AAQQES ETTQRAEREVTRMVIIMVIGFLVCWLPYASVAWFIFTHQGSEFGPLFMTIPAFF AKSSSI YNPMIYICMNKQFRHCMITTLFCGKNPFEGEEEGASSTKTEASSASSVSPA
Source	in vitro E.coli expression system
Target Names	rho
Protein Names	Recommended name: Rhodopsin
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.