



Recombinant Neoniphon sammara Rhodopsin (rho)

Product Code	CSB-CF019681NBW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79812
Product Type	Transmembrane Protein
Immunogen Species	Neoniphon sammara (Spotfin squirrelfish) (Holocentrus sammara)
Sequence	MNGTEGPYFYVPMVNTTGVVRSPEYYPQYYLVNPAFAVLGAYMFFLIIFGFPI NFLTLY VTLEHKKLRTPLNYILLNLAVADLFMVIGGFTTTMYSSMHGYFVLGRLGCNLEG FSATLG GMISLWSLAVLAIERWVVVCKPTSNFRFGENHAIMGVSLTWTMALACTVPPLV GWSRYIP EGMQCSCGIDYYTRAEGFNNEFVLYMFFCHFVPLIIFFCYGRLLCAVKEAA AAQQES ETTQRAEREVTRMVILMVIGYLVCPYASVAWFIFTHQGSEFGPLFMTIPAFF AKSSSI YNPVIYICMNKQFRNCMITTLFCGKNPFEGEEEGASSTKTEASSASSVSPA
Source	in vitro E.coli expression system
Target Names	rho
Protein Names	Recommended name: Rhodopsin
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.