



Recombinant Mouse Cysteinyl leukotriene receptor 1 (Cysltr1)

Product Code	CSB-CF006465MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q99JA4
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MYLQGTKQTFLENMNGTENLTTSLINNTCHDTIDEFRNQVYSTMYSVISVVGFF GNSFVL YVLIKTYHEKSAFQVYMINLAIADLLCVCTLPLRVVYYVHKGKWLFGDFLCRLTT YALYV NLYCSIFFMTAMSFRCVAIVFPVQINLVTQKKARFVCIGIWIFVILTSSPFLMY KSYQ DEKNNTKCFEPPQNNQAKKYVLILHYVSLFFGFIIPFVTIIVCYTMIILTLLKNTMK KNM PSRRKAIGMIIVVTA AFLV SFMPYHIQRTIHLHLLHSETRPCDSVLRMQKSVVITL SLAA SNCCFDPLLYFFSGGNFRRRLSTFRKHSLSMITYVPKKKASLPEKGEEICNE
Source	in vitro E.coli expression system
Target Names	Cysltr1
Protein Names	Recommended name: Cysteinyl leukotriene receptor 1 Short name= CysLTR1 Alternative name(s): Cysteinyl leukotriene D4 receptor Short name= LTD4 receptor
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.