



# Recombinant *Macaca fascicularis* C-X-C chemokine receptor type 4 (CXCR4)

<b>Product Code</b>	CSB-CF630683MOV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q28474
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Macaca fascicularis</i> (Crab-eating macaque) ( <i>Cynomolgus</i> monkey)
<b>Sequence</b>	MEGISIYTSNDNYTEEMGSGDYDSIKEPCFREENAHFNRIFLPTIYSIIFLTGIVGN GLVI LVMGYQKKLRSMTDKYRLHLSVADLLFVITLPFWAVDAVANWYFGNFLCKAVH VIYTVNL YSSVLILAFISLDRYLAIVHATNSQKPRKLLAEKVVYVGVWIPALLLTIPDFIFASV SEA DDRYICDRFYPNDLWVVVFQFQHIMVGLILPGIVILSCYCIISKLSHSHKQKQK ALKT TVILILAFFACWLPYYIGISIDSFILLEIHKQGCEFENTVHKWISITEALAFFHCCLNP I LYAFLGAKFKTSAQHALTSVSRGSSLKILSKGKRGHSSVSTESSESSSFHSS
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	CXCR4
<b>Protein Names</b>	Recommended name: C-X-C chemokine receptor type 4 Short name= CXC-R4 Short name= CXCR-4 Alternative name(s): Fusin Leukocyte-derived seven transmembrane domain receptor Short name= LESTR Stromal cell-derived factor 1 r
<b>Expression Region</b>	1-352
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.