



Recombinant Rabbit C-C chemokine receptor type 5 (CCR5)

Product Code	CSB-CF633610RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1ZY22
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	MDYSMSTALYDIDYGMSEPCQKIDVKQVAARLLPPLYSLVFIFGFVGNLLVVLIL ITCKK LKSMTDIYLLNLAISDLLFLLTLPLWAHYAAA EWDFGGAMCKVFTGMYHMGYF GGIFFII LLTIDRYLAIVHAVFALKARTVTFGVVTSVGTWVAAILVSLPDIIFTRSQKEGFRC SCSP HFPASQYQFWKNFHTIMRNILSLVPLLMIVCYSGILKTLLRCRNEKRRHRAV RLIFAI MVVYFLWAPYNVLLLNTFQEFFGLNNCSSSNRLDRAMQVTETLGMTHCCI NPVVYAFV GEKFRSYLSAFFRKHVAKRLCKHCPLLPRETPEPASSVYTRSTGEQEISVGL
Source	in vitro E.coli expression system
Target Names	CCR5
Protein Names	Recommended name: C-C chemokine receptor type 5 Short name= C-C CKR-5 Short name= CC-CKR-5 Short name= CCR-5 Alternative name(s): MIP-1 alpha receptor CD_antigen= CD195
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.