



Recombinant Chlorocebus aethiops C-X-C chemokine receptor type 4 (CXCR4)

Product Code	CSB-CF886927DSU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9TSQ8
Product Type	Transmembrane Protein
Immunogen Species	Chlorocebus aethiops (Green monkey) (Cercopithecus aethiops)
Sequence	<p>MEGISIYTSNDNYTEEMGSGDYDSIKEPCFREENAHFNRIFLPTIYSIIFLTGIVGN GLVI LVMGYQKKLRSM TDKYRLHLSVADLLFVITLPFWAVDAVANWYFGNFLCKAVH VIYTVNL YSSVLILAFISLD RYLAIVHATNSQKPRKLLAEKV VYVGVWIPALLLTIPDFIFASV SEA DDRYICDRFYPNDLWVVVFQFQHIMVGLILPGIVILSCYCIISKLSHSHKGHQK GK ALKT TVILILAFFACWLPYYIGISIDSFILLEI IKQGCEFENTVHKWISITEALAFFHCCLNP I LYAFLGAKFKTSAQH ALTSVSRGSSLKILSKGKRGGHSSVSTESESSSFHSS</p>
Source	in vitro E.coli expression system
Target Names	CXCR4
Protein Names	Recommended name: C-X-C chemokine receptor type 4 Short name= CXC-R4 Short name= CXCR-4 Alternative name(s): CD_antigen= CD184
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.