



# Recombinant Pongo abelii C-C chemokine receptor type 5 (CCR5)

<b>Product Code</b>	CSB-CF004844PYX
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P61756
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Sequence</b>	MDYQVSSPTYDIDYYTSEPCQKINVKQIAARLLPPLYSLVFIFGFVGNMLVILILIN CKR LKSMTDIYLLNLAISDLFFLLTVPFWAHYAAAQWDFGNTMCQLLTGLYFIGFFS GIFFII LLTIDRYLAIVHAVFALKARTVTFGVVTSVITWVVAVFASLPGIIFTRSQKEGLHY TCSS HFPYSQYQFWKNFQTLKIVILGLVLP LLVMVICYSGILKTL LRCRNEK KRHRAVR LIFTI MIVYFLFWAPYNI VLLLNTFQE FFLNCS S NR LDQAMQV TETLGMTHCCINP IYAFV GEKFRNYLLVFFQKHIAKRFCKCCSIFQQEAPERASSVYTRSTGEQEISVGL
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	CCR5
<b>Protein Names</b>	Recommended name: C-C chemokine receptor type 5 Short name= C-C CKR-5 Short name= CC-CKR-5 Short name= CCR-5 Short name= CCR5 Alternative name(s): CD_antigen= CD195
<b>Expression Region</b>	1-352
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.