



Recombinant *Nomascus leucogenys* C-C chemokine receptor type 5 (CCR5)

Product Code	CSB-CF004844NHM
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O97883
Product Type	Transmembrane Protein
Immunogen Species	<i>Nomascus leucogenys</i> (Northern white-cheeked gibbon) (<i>Hylobates leucogenys</i>)
Sequence	MDYQVSSPTYDIDYDTSEPCQKINVKQIAARLLPPLYSLVFIFGFVGNMLVILVLI NCKR LKSMTDIYLLNLAISDLFFLLTVPFWAHYAAAQWDFGNTMCQLLTGLYFIGFFS GIFFII LLTIDRYLAIVHAVFALKARTVTFGVVTSVITWVVAVFASLPGIIFTRSQKEGLHY TCSS HFPYSQYQFWKNFQTLKIVILGLVPLLVVICYSGILKTLRLCRNEKKRHRRAVR LIFTI MIVYFLFWAPYNIVLLLNTFQEFFGLNNCSSSNRLDQAMQVTETLGMTHCCINP IYAFV GEKFRNYLLVFFQKHIKHFCKCCSIFQQEAPERASSVYTRSTGEQEISVGL
Source	in vitro E.coli expression system
Target Names	CCR5
Protein Names	Recommended name: C-C chemokine receptor type 5 Short name= C-C CKR-5 Short name= CC-CKR-5 Short name= CCR-5 Short name= CCR5 Alternative name(s): CD_antigen= CD195
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.