



Recombinant Pygathrix bieti C-C chemokine receptor type 5 (CCR5)

| | |
|--------------------------|---|
| Product Code | CSB-CF004844FHN |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | O97880 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Rhinopithecus bieti (Black snub-nosed monkey) (Pygathrix bieti) |
| Sequence | MDYQVSSPTYDIDYYTSEPCQKVNKQIAARLLPPLYSLVFIFGFVGNILVVLILINCKR LKSMTDIYLLNLAISDLFFLLTVPFWAHYAAAQWDFGNTMCQLLTGLYFIGFFSGIFFII LLTIDRYLAIVHAVFALKARTVTFGVVTSVITWVVAVFASLPGIIFTRSQREGLHYTCSS HFPYSQYQFWKNFQTLKIVILGLVLP LLVMVICYSGILKTL LRCRNEK KRHRAVRLIFTI MIVYFLFWAPYNIVLLLNTFQEFFGLNNCSSSNRLDQAMQVTETLGMTHCCINPIIYAFV GEKFRNYLLVFFQKHI AKRFCKCCYIFQQEAPERASSVYTRSTGEQEISVGL |
| Source | in vitro E.coli expression system |
| Target Names | CCR5 |
| Protein Names | Recommended name: C-C chemokine receptor type 5 Short name= C-C CKR-5 Short name= CC-CKR-5 Short name= CCR-5 Short name= CCR5 Alternative name(s): CD_antigen= CD195 |
| Expression Region | 1-352 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |