



Recombinant Gorilla gorilla gorilla C-X-C chemokine receptor type 2 (CXCR2)

Product Code	CSB-CF636722GGZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q28422
Product Type	Transmembrane Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Sequence	<p>FNMESDSFEDFWKGEDLSNYSYSSALPPFLLDASPCEPESLEINKYFVVIYALV FLLSL LGNSLVILVILYSRVGRSVTDVYLLNLALADLLFALTLPIWAASKVNGWIFGTFLC KVVS LLKEVNFYSGILLACISVDRYLAIVHATRTLQKRYLVKFICLSIWGLSLLLALPV LLF RRTIYPSNVSPVCYEDMGNNTANWRMLLRILPQSFGFIVPLLIMLFCYGFTLRT LFKAHM GQKHRAMRVIFAVVLIFLLCWLPYNLVLLADTLMRTQVIQETCERRNHINQALD ATEILG ILHSCLNPLIYAFIQKFKCHGLLKILAIHGLISKDSLPKDSRPSFVGSSSGHT</p>
Source	in vitro E.coli expression system
Target Names	CXCR2
Protein Names	<p>Recommended name: C-X-C chemokine receptor type 2 Short name= CXC-R2 Short name= CXCR-2 Alternative name(s): High affinity interleukin-8 receptor B Short name= IL-8R B CD_antigen= CD182</p>
Expression Region	1-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>