



Recombinant Pig B1 bradykinin receptor (BDKRB1)

Product Code	CSB-CF816835PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8HZN9
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MASQTLVVFQASNQSQLPPPNATLCDGAQEAWHLLHKVLPTCVVAICSGGLL GNLFVLSV FLVPRRRLNAAEIYLAHLAASDLVFALGLPFWAETIRNGFHWPFGAPLCRVVN GVIKANL FISIFLVVAISRDRYRALVHPVASWRRRRRRHWAQATCVLIWTAGGLLSIPTFLL RSVQV VPELNVSACVLPFPHEAWAFVRTVELNVLGFLPLAAILFFNYHILAALRGREQL SRTRC GGPRDGKTTALILTLVAVFLLCWTPYHVCAFLEFLLHVRAIRGCFWEDFTDLGL QYTNFF AFINSCLNPVIYVFWGQLFRTKIWELYHRCLPRKLTAVSSSRKEIFQIFWRN
Source	in vitro E.coli expression system
Target Names	BDKRB1
Protein Names	Recommended name: B1 bradykinin receptor Short name= B1R Short name= BK-1 receptor
Expression Region	1-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.