




Recombinant Mouse C-C chemokine receptor type 8 (Ccr8)

Product Code	CSB-CF004847MO
Abbreviation	Recombinant Mouse Ccr8 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P56484
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose
Product Type	Transmembrane Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MDYTMENVTMTDYYPDFFTAPCDAEFLLRGSMLYLAILYCVLFVLGLLGNSLV ILVLVGCKKLRSITDIYLLNLAASDLLFVLSIPFQTHNLLDQWVFGTAMCKVVS G LYYIGFFSSMFFITLMSVDRLAIVHAVYAIKVRTASVGTALSLTVWLA AVTATIP LMVFYQVASEDGMLQCFQFYEEQSLRWKLFTHFEINALGLLLPFAILLFCYVRIL QQLRGCLNHNRTRAIKLVLTVVIVSLLFWVPFNVALFLTSLHDLHILDGCATRQR LALAIHVTEVISFTHCCVNPVIYAFIGEKFKKHLMDFVQKSCSHIFLYLGRQMPV GALERQLSSNQRRSSHSSTLDDIL
Research Area	Immunology
Source	in vitro E.coli expression system
Target Names	Ccr8
Protein Names	Recommended name: C-C chemokine receptor type 8 Short name= C-C CKR-8 Short name= CC-CKR-8 Short name= CCR-8 Alternative name(s): CD_antigen= CDw198
Expression Region	1-353aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	42.8 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.