



Recombinant Danio rerio Membrane progesterin receptor alpha-B (paqr7b)

Product Code	CSB-CF801898DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q801G2
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MATVVMEQIGRLFINAQQLRQIPRFLESAFPKLPTVMVSDVPWVFRESHIITG YRPPDQ NWRYYFLTLFQRHNESVNVWTHLLASLIILVKFQELSETVDFLRDPHAQPMFIL LLAAFT YLGCSALAHLLSAKSEISHYTFYFLDYVGVAVYQYGSALAHFYVVEEEWHAQ VRTFFLP ASAFLAWLSCTGCCYGKYASPCLKFVHKLFQVPSGLAYCLDISPVLHRIYRC YSSEHW CADQAVVYHCYQVLFFLISAYFFSYPHPERWFPGRCDFIGQGHQIFHVFLVLC TLVQIEA VRLDYTERRRLYEHLHGDLAHDVAVALFIFTACCSALTAFYVRKRKVKTYLEEKQE
Source	in vitro E.coli expression system
Target Names	paqr7b
Protein Names	Recommended name: Membrane progesterin receptor alpha-B Short name= mPR alpha Alternative name(s): Progesterin and adipoQ receptor family member VII, b
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.