



Recombinant Rana pipiens Rhodopsin (RHO)

Product Code	CSB-CF019681RJN
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31355
Product Type	Transmembrane Protein
Immunogen Species	Lithobates pipiens (Northern leopard frog) (Rana pipiens)
Sequence	MNGTEGPNFYIPMSNKTGVVRSPFDYPQYYLAEPWKYSVLAAYMFLILLGLPI NFMTLY VTIQHKKLRTPLNYILLNLGVCNHFMLVLCGFTITMYTSLHGYFVFGQTGCYFEG FFATLG GEIALWLSLVVLAIERIVVCKPMSNFRFGENHAMMGVAFTWIMALACAVPPLF GWSRYIP EGMQCSCGVDYITLKPEVNNESFVIYMFVVHFLIPLIIISFCYGRLVCTVKEAAA QQQES ATTQKAEKEVTRMVIIMVIFFLICWVPYAYVAFYIFTHQGSEFGPIFMTVPAFFAK SSAI YNPVIYIMLNKQFRNCMITTLCCGKNPFGDDDASSAATSKTEATSVSTSQVSPA
Source	in vitro E.coli expression system
Target Names	RHO
Protein Names	Recommended name: Rhodopsin
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.