



Recombinant *Macaca fascicularis* G-protein coupled receptor 1 (GPR1)

Product Code	CSB-CF856817MOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q95LH1
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Sequence	MEDLEETLFEEFENYSYALDYYSLESDLEEKVQLGVVHVVSLVLYCLSFVLGIP GNAIVI WFTGFKWKRTVSTLWFLNLAIADFIFLLFLPLYISYVVMNFHWPFGIWLCKANS FTAQLN MFASVFFLTVISLDHYIHLIHPVLSHRHRTLKNSLIVIFIWLLASLIGGPALYFRDT VE FNNHTLCYNNFQKHDPDLTVIRHHVLTWVKYIVGYLFPLLTMSICYLCLILKVKK RSILI SSRHFWTILAVVVAFVVCWTPYHLFSIWELTIHHNSYSHHVMQAGIPLSTGLAF LNSCLN PILYVLISKKFQARFRSSVAEILKYTLWEVSCSGTVSEQLRNSETKNLCLLETAQ
Source	in vitro E.coli expression system
Target Names	GPR1
Protein Names	Recommended name: G-protein coupled receptor 1
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.