



# Recombinant Human CX3C chemokine receptor 1 (CX3CR1)

<b>Product Code</b>	CSB-CF006236HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P49238
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	MDQFPESVTENFEYDDLAEACYIGDIVVFGTVFLSIFYSVIFAIGLVGNLLVVFAL TNSK KPKSVTDIYLLNLALSDLLFVATLPFWTHYLINEKGLHNAMCKFTTAFFFIGFFG SIFFI TVISIDRYLAIVLAANSMNRTVQHGVTVISLGVWAAAILVAAPQFMFTKQKENE CLGDYP EVLQEIWVPLRVNVTNFLGFLPLLIMSICYFRIIQTLFSCKNHKKAKAIKLILLV IVF FLFWTPYNYMIFLETCLKLYDFFPSCDMRKDLRLALSVTETVAFSHCCLNPLIYAF AGEKF RRYLYHLYGKCLAVLCGRSVHVDFSSSESQRSRHSVLSNFTYHTSDGDAL LLL
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	CX3CR1
<b>Protein Names</b>	Recommended name: CX3C chemokine receptor 1 Short name= C-X3-C CKR-1 Short name= CX3CR1 Alternative name(s): Beta chemokine receptor-like 1 CMK-BRL-1 Short name= CMK-BRL1 Fractalkine receptor G-protein coupled rece
<b>Expression Region</b>	1-355
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.