



Recombinant Drosophila melanogaster Protein rhomboid (rho)

Product Code	CSB-CF323222DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P20350
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MENLTQNVNETKVDLGQEKEKEASQEEEHATAAKETIIDIPAACSSSSNSSSYD TDCSTA SSTCCTRQGEHIYMQREAIPATPLPESEDIGLLKYVHRQHWPWFILVISIIEIAIFA YDR YTMPAQNFGLPVPIPSDSVLVYRPDRRLQVWRFFSYMFLHANWFHLGFNIVIQ LFFGIPL EVMHGTARIGVIYMAGVFAGSLGTSVVDSEVFLVGASGGVYALLAAHLANITLN YAHMKS ASTQLGSVVIFVSCDLGYALYTQYFDGSAFAKGPQVSYIAHLTGALAGLTIGFL VLKNFG HREYEQLIWWLALGVYCAFTVFAIVFNLINTVTAQLMEEQGEVITQHLLHDLGV S
Source	in vitro E.coli expression system
Target Names	rho
Protein Names	Recommended name: Protein rhomboid EC= 3.4.21.105 Alternative name(s): Protein veinlet
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.