



Recombinant Rabbit C-X-C chemokine receptor type 1 (CXCR1)

Product Code	CSB-CF325015RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P21109
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	MEVNVWNMTDLWTFEDEFANATGMPPVEKDYSPCLVVTQTLNKYVVVVIYA LVFLLSLL GNSLVMLVILYSRSNRSVTDVYLLNLAMADLLFALTMPIWAVSKEKGWIFGTPL CKVVSL VKEVNFYSGILLACISVDRYLAIVHATRRTLQKRHLVKFICLGIWALSILSLPFF LFR QVFSPNNSPVCYEDLGHNTAKWRMVLRLPHTFGFILPLLVMLFCYGFRTLRTL FQAHMG QKHRAMRVIFAVVLIFLLCWLPYNLVLLADTLMRTHVIQETCQRRNDIDRALDA TEILGF LHSCLNPIIYAFIQNFRNGFLKMLAARGLISKEFLTRHRVTSYTSSSTNVPSNL
Source	in vitro E.coli expression system
Target Names	CXCR1
Protein Names	Recommended name: C-X-C chemokine receptor type 1 Short name= CXC-R1 Short name= CXCR-1 Alternative name(s): High affinity interleukin-8 receptor A Short name= IL-8R A CD_antigen= CD181
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.