



Recombinant Rabbit CX3C chemokine receptor 1 (CX3CR1)

Product Code	CSB-CF644447RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2KTE1
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	MTTLYSDWATESFEYDESSEACFIGDIVAFGTIFLSIFYSLVFAFGLVGNLLVVC ALTSS RKPKSITDIYLLNLALS DLLFVATLPFWTHYVISEQGFHNAVCKLTTALFFIGFFG GIFF ITVISIDRYMAIVLAANSINNRTVQHGVTTSLGVWAAAILVAAPQFMFTKQKGNE CLGDY PEVLQDIWPVLRNTEANFLGFLLPVLIMSYCYFRIIQLFSCKNHKKAKAIKLILL VIV FFLFWTPYNVMIFLETCLKLYGFFPNCMDMKRDLRLALSVTETVAFSHCCLNPLIY AFAGQK FRRYLRHLSRKQCQAVLCGRP VHVSFSPSESQRSRQESIVSSNFTHYTS DGDA SLLL
Source	in vitro E.coli expression system
Target Names	CX3CR1
Protein Names	Recommended name: CX3C chemokine receptor 1 Short name= C-X3-C CKR-1 Short name= CX3CR1 Alternative name(s): Fractalkine receptor
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.