



# Recombinant *Aedes aegypti* Elongation of very long chain fatty acids protein AAEL008004 (AAEL008004)

<b>Product Code</b>	CSB-CF624917AXQ
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q1HRV8
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Aedes aegypti</i> (Yellowfever mosquito) ( <i>Culex aegypti</i> )
<b>Sequence</b>	MALIMKYIDSMHHYMDKYGDPRTKDWPLMSSPFPTLALCLGYVYLVKVLGPRL MENRKPF QLRNTLILYNFVQVVFSAWLFYEIGISGWLTGHYNFRCQPVDYSNHPKTLRMV HACWWYY FSKFTEFFDTFFFVMRKKTSQVSTLHVIIHHGCMPSVWFGVKFTPGGHSTFF GLLNTFVH IVMYTYYLFTAMGPQFQKYLWWKKYLTSLQMVQFVAIMVHAFQLLFIDCNYPK AFVWWIG MHAVMFLFLFNEFYQSTYKATKRRRAAAAEARRLAAEEAKLQNGSAVSSNGS AITANGHH GKNGSVHHHSNGSATSNGTSLLSNGVGSNKAADYYVRGDLPAEIEITQRQPS SRNQVQ
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	AAEL008004
<b>Protein Names</b>	Recommended name: Elongation of very long chain fatty acids protein AAEL008004 EC= 2.3.1.n8 Alternative name(s): 3-keto acyl-CoA synthase AAEL008004
<b>Expression Region</b>	1-358
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.