



Recombinant Serpentine receptor class epsilon-33 (sre-33)

Product Code	CSB-CF515719CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O18175
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	<p>MIINSNSSTIFSSIWLPVFFYVEPLDQQVIISILELMIYLVCIHLVNVSLHVALKIRLF H RNLYILALPMFGMWYELIIGKFITIAAYRLKLLGLDFELGEHTAIWTNDPGKVLLVA SLNG LELLIFGGFLQWHYMYSWIFGVLTVAVERVIASVLIENYESNTQNLMPAILLIISQ FLSI SMAFGLLFQKVGPLSAHFPWMISCPISVAAYVFVKKNESFRREIKNPGRKRIF TLSQQF QVKENLRVLHLGTRLVFAVLSFIGICGCGIAALHYKIVPSYYCHLIENVLFLNPFLI GLT AMLSIPQWKEQFMKSFLT VRLFRNRRKPVHIVVEIEECAKKKNDVETNLYFKQL ANSWI</p>
Source	in vitro E.coli expression system
Target Names	sre-33
Protein Names	Recommended name: Serpentine receptor class epsilon-33 Short name= Protein sre-33
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.