



Recombinant Dog C-C chemokine receptor type 3 (CCR3)

Product Code	CSB-CF714527DO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q64H34
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sequence	MEATTAEIKTMDESIQTTVFDYENSLPCEKVNKHLGAQFLPPLYSLVFIGLLG NWWW VILTKYKRLWIMTNIFLLNLAISDLLFLFTLVFWIHYTGWWDWVFGSRMCKLISGL YYLG LYGEIFFIILLTIDRYLAIVHAVFALRARTVTFGIITSVLTWGLAGLAALPEFIFHES QK ESEQFVCTPLYPKDQEDNWKRFBALRMNILGLALPLLIMVVCYSGIIKTLLRCP SKKKYK AIRLIFVIMVFFIFWTPYNLVLLSAFQTIFFETTCEQSKQLDVAMQITEVIAYTH CCI NPIIYAFVGERFRKHLCHFFRRNVATYLGKYIPFLPSEKLERSSSVSPSTGEQE FSFVF
Source	in vitro E.coli expression system
Target Names	CCR3
Protein Names	Recommended name: C-C chemokine receptor type 3 Short name= C-C CKR-3 Short name= CC-CKR-3 Short name= CCR-3 Short name= CCR3 Short name= CKR3 Alternative name(s): CD_antigen= CD193
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.