



Recombinant Human Olfactory receptor 5T2 (OR5T2)

Product Code	CSB-CF850889HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGG2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSYSIYKSTVNIPLSHGVVHSFCHNMNCNFMHIFKFLDFNMKNVTEVTLFVLK GFTDNL ELQTIFFFLFLAIYLF TLMGNLGLILVVIRDSQLHKPMYYFLSMLSSVDACYSSVI TPNM LVDFTTKNKVISFLGCVAQVFLACSFSGTTECFLLAAMAYDRYVAIYNPLLYSVS MSPRVY MPLINASYVAGILHATIHTVATFSLFCGANEIRRVFCDIPPLLAISYSDTHTNQL LLFY FVGSIELVTILIVLISYGLILLAILKMYSAEGRRKVFSTCGAHLTGVSIIYG TILFMY VR PSSSYASDHDMIVSIFYTIVIPLLNPVIYSLRNKDKDSMCKMFGKNQVINKVYF HTKK
Source	in vitro E.coli expression system
Target Names	OR5T2
Protein Names	Recommended name: Olfactory receptor 5T2 Alternative name(s): Olfactory receptor OR11-177
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.