



Recombinant *Macaca mulatta* (Rhesus macaque) G-protein coupled receptor 15 (GPR15)

Product Code	CSB-CF009764MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O97663
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca mulatta</i> (Rhesus macaque)
Sequence	MDPEETSVYLDYYYATSPNPDIRETHSHVPYTSVFLPVFYTAVFLTGVLGNLVL MGALHF KPGSRRLIDIFIINLAASDFIFLVTLPLWVDKEASLGLWRTGSFLCKGSSYMISVN MHCS VFLLCMSVDRYLAIVCPVVS RKFRRTDCAYVVCASIWFISCLLGLPTLLSRELT LIDDK PYCAEKKATPLKLIWLSLVALIFTFFVPLLSIVTCYCCIAARKLCAHYQQSGKHNNK LKKSI KIIFIVVAAFLVSWLPFNTFKLLAIVSGLQQERYFPSAMLQLGMEVSGPLAFANS CVNPF IYYIFDSYIRRAIVHCLCPCLKNYDFGSSTETSDSHLTKALSTFIHAEDFTRRRKR SVSL
Source	in vitro E.coli expression system
Target Names	GPR15
Protein Names	Recommended name: G-protein coupled receptor 15
Expression Region	1-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.