



# Recombinant Bovine Proteinase-activated receptor 2 (F2RL1)

<b>Product Code</b>	CSB-CF643165BO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q2HJA4
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Sequence</b>	SLIGNVDNSPVVAGRGVTVKPGFSVDEFSTSVLTGKLTTVFLPVVYTIVFVVG LPSNGMA LWVFLFRTKKKHPAVIYMANLALADLLSVTWFPKIAHYHIHGNNWIYGESLCKVL IGFFY GNMYCSILFMTCLSVQRYWVIVNPMVHPKKQANIAIGVSLGIWLLILLTIPLYV VYKQTS YIRALNITTCHDVLPEEVLVGDMFNYFLSLAIGVFLFPAFLTASAYVLMIRTLQSS AMDE SSGKKRRRAIKLIVTVLAMYLICFTPSNLLLHVHYFLIKTRGQSHVYALYIVALCL STLN SCIDPFVYYFISQDFRDHAKNALLCRSVRTVKRMQVSLSSKKFSGKSSSYSS STS VKGS Y
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	F2RL1
<b>Protein Names</b>	Recommended name: Proteinase-activated receptor 2 Short name= PAR-2 Alternative name(s): Coagulation factor II receptor-like 1 Thrombin receptor-like 1
<b>Expression Region</b>	35-395
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.