



Recombinant Mouse P2Y purinoceptor 4 (P2ry4)

Product Code	CSB-CF888329MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9JJS7
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MTSADSLIFTSLGPSPSSGDGCKFNEEFKFILLPLSYAVVFLGLALNAPTLW LFLFRL RPWDATATYMFHLALSDTLYVLSLPTLVYVYAARNHWPFGTGFCKFVRFIFY WNLYCSVL FLTCISVHRYMGICHPLRAIRWGRPRFAGLLCLGVWLVVAGCLVPNLFFVTTNA NGTTIL CHDTTLPEEFDHYVYFSSTIMVLLFGFPFLITLVCYGLMARRLYRPLPGAGQSS SRLRSL RTIAVVLTVFAVCFVPFHITRTIYLLARLLNAECRVLNIVNVVYKVTRPLASANSC LDPV LYLFTGDKYRNQLQQLCRGSTPKRRTTASSLALVTLHEESISRWADIHQDSIFP AYEGDR L
Source	in vitro E.coli expression system
Target Names	P2ry4
Protein Names	Recommended name: P2Y purinoceptor 4 Short name= P2Y4
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.