



Recombinant Mouse Probable G-protein coupled receptor 52 (Gpr52)

Product Code	CSB-CF009818MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C5J4
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MNESRWTEWRILNMSSSIVNVSEHHSCPLGFGHYSVEDVCIFETVVIVLLTFLII SGNLT VIVFVHCAPLLHHYTTTSYFIQTMAYADLLVGVTCCLVPTLSLLHYSTGVHESLTCQ VFGYI ISVLKSVSMACCLACISVDRYLAITKPLSYNQLVTPCRLRICIIMIWIYSCLIFLPSFF GW GKPGYHGDIFEWCATSWLTSAYFTCFIVCLLYAPAALVVCFTYFHIFKICRQHT KEINDR RARFPSHEVEASREAGHSPDRRYAMVLFRITSVFYMLWLPYIIYFLESSRVLD NPTLSF LTTWLAISNSFCNCVIYSLSNSVFRLGLRRLSETMCTSCVCAKDQEAQDPKPR RRANSCS I
Source	in vitro E.coli expression system
Target Names	Gpr52
Protein Names	Recommended name: Probable G-protein coupled receptor 52
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.