



Recombinant Human G-protein coupled receptor 183 (GPR183)

Product Code	CSB-CF009792HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P32249
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDIQMANNFTPPSATPQGNDLCDLYAAHSTARIVMPLHYSLVFIIGLVGNLLALV VIVQNR KKINSTTLYSTNLVISDILFTTALPTRIAYYAMGFDWRIGDALCRITALVFYINTYA GVN FMTCLSIDRFIAVVHPLRYNKIKRIEHAKGVCIFVWILVFAQTLPLLINPMSKQEA ERIT CMEYPNFEETKSLPWILLGACFIGYVLPLIILICYSQICCKLFRITAKQNPLTEKSG VNK KALNTIILIVVFLCFTPYHVAVIIQHMIKKLRFNSNFLECSQRHSFQISLHFTVCLM NFN CCMDPFIYFFACKGYKRKVMRMLKRQVSVSISAVKSAPEENSREMTETQMM IHSKSSNG K
Source	in vitro E.coli expression system
Target Names	GPR183
Protein Names	Recommended name: G-protein coupled receptor 183 Alternative name(s): Epstein-Barr virus-induced G-protein coupled receptor 2 Short name= EB12 Short name= EBV-induced G-protein coupled receptor 2
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.