



Recombinant Serpentine receptor class X 45 (srx-45)

Product Code	CSB-CF330360CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34490
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	<p>MFQILMENVEVQHIDRIAALLIFFTSFIGFACNTFIAFYIRRLSLLRNSFGRLQLQ AAG DAVFLVWAFYFAPVLFFDIKPLQSLAIAARFAQLCLICYDISIYTHLVISLNRFISL YF PTSYQNIFTERFTTFLICSIIFVSGFSWFLVIRDCQMGFSIPRWMLDYVSPPE MINVY YAEFFRGLIVISMFAITNSFTFCRMHMHNRKKQTATVFETTQQKKRAVETRFV QQVTMQ GLLYVIELVTYFYISLRFPVPLEPVELAKSSNRWPNFLLTTYAWILVHALDGVITLI FNK QFRSVLRHPCRSQEALNISKTPSRRSRFTTNRDESNRKKSSILACTFISNSNTG YGSSVH V</p>
Source	in vitro E.coli expression system
Target Names	srx-45
Protein Names	Recommended name: Serpentine receptor class X 45 Short name= Protein srx-45
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.